

TECHNICAL SPECIFICATIONS

INDICATION	DVT prophylaxis, reduction of post-operative pain and swelling, reduction of wound healing time
SIZE	5.3 x 5.3 x 2.4 in (135 x 135 x 60 mm)
WEIGHT	730 g / 1.6 lb (including battery)
POWER INPUT	9.7 VDC
BATTERY*	7.2 V (rechargeable NiMH)
Battery life	5h (after full recharge)
Recharge time	4-5 h
CHARGER	Built-in
MAINS SUPPLY	100-240V, 50/60 Hz
INFLATION MODE	Sequential compression
PRESSURE (preset)	
Calf	50 mmHg (± 10%)
Foot	130 mmHg (± 10%)
Thigh	50 mmHg (± 10%)
COMPRESSION	
DURATION	
Calf	Up to 10 sec
Foot	Up to 6 sec (depending on model)
Thigh	Up to 15 sec (depending on model)
PULSE FREQUENCY	
Calf	1 per minute
Foot	2 per minute
Thigh	1 per minute

*Battery life and recharge may vary depending on usage.

WARNING: Use only the battery pack provided by MCS.

Use of any alternative power source could cause irreparable damage to the system!

All MCS products are LATEX FREE
FDA CLASS II Prescription device